

SURVEY REPORT

RS.....

INDONESIAN NATIONAL DIAGNOSTIC REFERENCE LEVEL SURVEY

HEAD INFO			
PRACTICE NAME	RS.....	AGE OF GROUP
PROTOCOL	SCANNER

SETTINGS			
kVp	Rotation time (s)
mAs	No. of Phases	1; 2; 3; 4; 5; 6;
Pitch	Helical/Axial	Helical / Axial
Contrast	Yes / No	Detector Configuration x (.....)
Dose Modulation	Yes / No	Iterative Constructive	Yes / No
CTDI	CTDIvol / CTDIw		

Comments

SURVEY DATA					
Patients	Average CTDI (mGy)	Total DLP (mGy.cm)	Weight (Kg)	Age (years)	Sex (F/M)
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					

Ket:

- No of Phase tidak termasuk scano/scout/topo
- Average CTDI tidak termasuk scano/scout/topo, = (CTDI 1 + CTDI 2 + CTDI 3 +)/No of Phase
- Nama Protocol = 1. Head, 2. Chest, 3. Neck, 4. Abdopelvis, 5. CestAbdoPelvis, 6. Lumbar Spine
- Age Group = 1. Baby (0-4 years), 2. Children (5-14years), 3. Adult (>15 years)
- nilai mAs dapat diisi dengan mA setting (mA awal) jika menggunakan dose modulation.